

Powerful filtering, excellent permeability

AirQUEEN Nano Mask

Masks are a must!

Introducing easy-to-breathe mask with exceptional filter





MADE IN
K  **REAA**

Made with our own proud
technology-AirQUEEN Nano Mask

- Nanofiber Filter Manufactured by: LEMON Co., Ltd.
- Mask Manufacture & Distribution: TOPTEC Co., Ltd.

Breathing New Technology, Nanofiber Nanofiber Filter

Nano fibers are only 1 nanometer, 1 billionth of a meter, in thickness
Nanofiber is a new material made by building these fibers sterically
for a fishnet structure that has better permeability, blockability, and
durability than any fiber ever invented.



Permeability



Blockability




Durability

Nanofiber Filter has a different fiber thickness / fiber thickness / filtering method from the existing electrostatic filter

Melt-Blown Filter

5,000 times enlarged


The fiber is thick
Average space between the
fibers is over 10 μ m



Nanofiber Filter

15,000 times enlarged

Very thin, nano thickness
Average space between the fibers
is 1.0 μ m, with self-filtering




The Nanofiber Filter from is a cutting-edge technology applied to global brands and products

It is exclusively distributed to **THE NORTH FACE** as the brand **FUTURELIGHT**.

It is also used for **RENERGIE mask pack** from **LANCÔME** that touches the skin directly.

This technology is used for **smart devices** that sends **sounds** but needs **protection against dust and moisture**.



We asked professor Kim Iksu, the most renowned professor in the fiber research area.

What kind of material is **Nanofiber**?

The thickness of a nano fiber is only a thousandth of a strand of hair.
Nanofiber is a material made by building these fibers sterically for a fishnet structure,
and nanofiber of 85% processing rate has both permeability and blockability.

The Nanofibers of LEMON have consistent thickness and pore size.
This means the ratio of the pores in the general area becomes very high,
which is relevant to permeability.

Because it touches the **nose** and the **mouth**,
this fiber is **very appropriate for masks**
where permeability is important

[From interview with professor Kim Iksu of Department of Textile at Shinshu University, Japan]





**Long-lasting freshness
from morning to evening!**

The effectiveness of electrostatic masks decreases with each breath due to the moisture
Nanofiber Filter has great permeability, so you can feel fresh even after a long term use



Excellent filtering and permeability of Nanofiber

Nanofiber smart filter blocks pollutants and enables easier breathing

**If you care about the skin,
permeability is even more important.**

Wearing mask all day causes skin trouble (C Newspaper)

Must wear permeable masks (M Newspaper)

Pimple breaks out due to wearing masks
(National Health Insurance Service blog)

There is moisture in our respiratory system
If a mask has bad permeability, the moisture from within
cannot escape, which causes eczema

If you wear a mask for a long time,
you really need a permeable mask

Because we need to wear it all day,
we made it lighter.

AirQUEEN Nano mask 3.71g


Lighter than a sheet of paper, so you won't even realize wearing it!



※ a sheet of paper = 4.72g

3D ergonomic design to fit perfectly on all face types.

Nose bridge prevents sliding or fogging up.





**Individually wrapped
for easy carrying and using**





(19) 대한민국특허청(KR)

(12) 등록특허공보(BI)

(45) 공고일자 2011년06월09일

(11) 등록번호 10-1040055

(24) 등록일자 2011년06월02일

(51) Int. Cl.
D01D 5/00 (2006.01) **D01D 13/00** (2006.01)
B82Y 40/00 (2011.01)

(21) 출원번호 10-2011-0016876

(22) 출원일자 2011년02월24일

심사청구일자 2011년02월24일

(30) 우선권주장 JP-P-2010-272071 2010년12월06일 일본(JP)

(56) 선행기술조사문헌 JP2009052163 A KR1020070026744 A KR1020070047282 A JP2008519175 A

전체 청구항 수 : 총 14 항 심사관 : 최봉문

(73) 특허권자 신승 다이가부

일본 나가노켄 마쓰모토시 아사히 3초메 1번 1교

주식회사 토크

경상북도 구미시 산동면 봉산리 366번지

(72) 발명자 이재관

경상북도 구미시 산동면 봉산리 366 (주)토크

김익수

T.3898567나가노현우에다시토키타3-15-1국립대학

법인신유대학신유학부소속

(74) 대리인 오광환


(54) 전계방사장치 및 나노섬유 제조장치

(57) 요약

노즐 블록을 집지 한 상태에서 컬렉터에 고전압을 인가하여 전계 방사를 하는 경우라도, 균일한 품질을 가지는 나노 섬유를 높은 생산성으로 대량생산 하는 것이 가능한 전계방사장치를 제공한다.

컬렉터(150)와, 노즐 블록(110)과, 정전극이 컬렉터(150)에 접속되고, 부전극이 노즐 블록(110)에 접속되는 동시에 해당 부전극의 전위가 집지 전위로 떨어진 장치(160)를 구비하는 전계방사장치(20)로서, 컬렉터(150)를 둘러싸는 위치에 회전 자유롭게 배치된 절연성 또한 다공성의 엔들리스 벨트로 이루어지는 보조 벨트(172) 및 해당 보조 벨트(172)를 강직시트(170)의 이송 속도에 대응하는 회전 속도로 회전시키는 보조 벨트 구동장치(174)를 가지는 보조 벨트 장치(170)를 추가로 구비하는 전계방사장치.

대표도



The Nanofiber Filter of AirQUEEN Nano Mask is manufactured with patent technology.


- Cutting-edge technology that mass-produces consistent nano fiber filters
- Selected as one of 100 technologies to lead Korea in 2020
- Approximately 50 patents around the world

- Patent Registration no. : 10-1040055
- Name of Invention: An electrospinning device and an apparatus for manufacturing nano-fiber

Nanofiber Filter with excellent permeability and perfect protection from water and dust, used in AirQUEEN Masks, is the fruit of LEMON, a company listed through KOSDAQ, researching and developing nano fiber technology for 14 years

The Nanofiber Filter of AirQUEEN Mask is produced in a clean and sanitary facility.

- Proven suitable for ISO 9001 quality management system
- Proven suitable for ISO 14001 environment management system



The Nanofiber Filter used for AirQUEEN Masks is produced in a facility approved for environment management system, so it is safe and sanitary.



Test Report No. F69101/LF-CTSAYGA18-00327 Issued Date: 2018. 01. 16. Page 1 of 17

Korea



The following sample(s) was/were submitted and identified by/on behalf of the client as:

SGS File No. : AVGA18-00327

Product Name : Nano membrane (sanitary pad)

Item/Part Name : N/A

Received Date : 2018. 01. 09

Test Period : 2018. 01. 09. ~ 2018. 01. 16.

Test Requested : One hundred-Seventy four (174) substances in the Candidate List of Substances of Very High Concern (SVHC) for authorization published by European Chemicals Agency (ECHA) on July 7, 2017 regarding Regulation (EC) No 1907/2006 concerning the REACH.
Eight (8) substances in the Public Consultation List of potential Substances of Very High Concern (SVHC) published by European Chemicals Agency (ECHA) on Sep. 5, 2017 regarding Regulation (EC) No 1907/2006 concerning the REACH.

Test Method : Please refer to next page(s).

Test Result(s) : Please refer to next page(s).

Summary : According to the specified scope and evaluation screening, the test results of SVHC are 0.1% (w/w) in the articles of the submitted sample.

SGS Korea Co., Ltd

Jeff Jang

Jeff Jang / Chemical Lab Mgr

This document is prepared by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at www.sgs.com/sgs/termsandconditions. This document is the property of SGS Korea Co., Ltd. and is not to be distributed outside the Company's premises without the prior written consent of SGS Korea Co., Ltd. The Company's sole responsibility is to its Client and this document does not constitute a contract. The Client is responsible for their rights and obligations under the transaction documents. This document cannot be reproduced, stored in a retrieval system or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, or by any information storage and retrieval system, without the prior written approval of SGS Korea Co., Ltd. The test results shown in this test report refer only to the samples.



SGS Korea Co., Ltd

300, The G Center, 78, Ulsan, Daejeon, Ansan, Gyeonggi-do, Korea 14117
T +82 213 4038 2000 F +82 2131 4039 2100 www.sgs.com/korea

Member of the SGS Group (Societe Generale de Surveillance)

The Nanofiber Filter of AirQUEEN Mask is approved by European SGS, REACH, and RoHS, and therefore is safe.

- Cadmium N/D
- Lead N/D
- Mercury N/D
- Hexavalent chromium N/D
- Other 182 categories N/D



Metal material and harmful compound not detected



Organic solvent and compound not detected

Organic solvent and compound not detected

Proven that metal material and harmful compound were not detected.



AirQUEEN Nano Mask is an industrial product that passed the tests of KEMTI and is [currently in process] of being approved by the MFDS


* Korea Environment Merchandise Testing Institute, KEMTI

AirQUEEN Nano Mask

is a product that applied for the KF approval at the MFDS early April of 2022 after passing the test done by KEMTI, a testing institute designated by the MFDS.

It is manufactured with the same technology as Technoweb Disposable dustproof mask (KF94) from FTENE, the subsidiary company of TOPTEC, but if the manufacturer or brand is changed, the approval from the MFDS needs to be renewed, and until then, it cannot be publically distributed.

However, after the KF certification, 80% of the supply will be distributed as public, so there might be a shortage



We are currently also waiting for the result of N95 and CE certification.

Quality Indication According to Electric Products and Household Safety Law

- Product: AirQUEEN Mask
- Manufactured in: Korea
- Manufacturer: TOPTEC
- Distributor: TOPTEC
- Weight: 4.38g
- Product Type: Mask
- Expiry Date: 36 months
- Storage Directions: closed container, in room temperature (1~30°C)
- Entire Ingredients:
felt (outer fabric, filter, inner fabric), plastic sheath wire,
polypropylene ring, nylon strap









AirQUEEN Mask is an industrial product.

AirQUEEN Mask uses the Nanofiber Filter technology of LEMON and is manufactured by TOPTEC.

Currently, the MFDS approval is being processed, and when it is approved as sanitation mask, 80% of the product supply will be distributed as public distribution.




Air QUEEN

The Vital Shield for You and Your Loved Ones



NANO FILTER TECHNOLOGY



Introducing Air QUEEN mask from Korea

It's a mask with a difference

The upcoming Korean brand that has manufactures its products under the guidance of researchers and analyses and understands the problems of its consumers. Air QUEEN tries to bring in the best of technology and information to provide wholesome benefit to customers.



Thin and light



Nano fiber filter



Ergonomic design

Experience the obvious difference in the queen mask


It covers from the nose to the chin **completely**, and efficiently blocks the inhaling of the fine-dust through the adhesion structure, also shows **high stability** when wearing one.

Completely blocks outside harmful substance through **Nano fabric filter**.


Comfortable wearing feeling by covering down to the chin through **3 layer structure three dimensional design**

Lighter than 1 piece of A4 paper!


MASK COMPARISON TEST




THE DIFFERENCE BETWEEN NANO FILTER AND NORMAL FILTER



Magnified x500 through microscope



Magnified x500 through microscope



THE DIFFERENCE BETWEEN NANO FILTER AND NORMAL FILTER

The pre-existing static filter

The pre-existing static filter uses the method where the dust sticks to the fabric due to static.

As a characteristic of the static, if moisture occurs due to breathing or is brought in from the outside, the filtering efficiency drops drastically.

Our company's Nano fiber filter

Fine dust has been filtered on the densely tangled Air Queen Nano fabric. It filters through a Nano fabric structured not instead of static, it maintains high filtering efficiency for 24 hours even when it contacts with moisture.

LIGHTER THAN ANY OTHER MASK

It is **light and comfortable all day**, as if you are not wearing one.



Pre-existing mask over 5g




Air queen mask, lighter than 1
piece of A4 paper! 4.38g
(A4 paper _ 4.72g)

The pre-existing Mask is **hard to breathe** in because it's too thick.

Next generation new material **Nano Membrane** that lets the air through

BREATHING FIBER




Most existing masks are less efficient due to moisture generated when breathing in an electrostatic manner.


NANO-membrane is 0.1 μm larger than air particles and smaller than harmful particles, and is filtered by the nanofiber structure itself, allowing for long wear and breathability, blocking, and durability.

‘Nano Membrane’

NICOTINE FILTERING TEST



Fill the smoke on the left cube after place the mask between two cubes



Vacume the smoke from the other side of cube





Able to see the nicotine on the nano filter



Face side is still clean after filtered from nanofilter

AIR VENTILATION TEST





STEP 1


Hang the ear band to ears and adjust dotted line to your jawline

STEP 2

Check if both jawline stick to your mask

STEP 3

Press the nose wire and make sure if it fits the face



Cautions

- Do not use in environment of oxygen concentration that is less than 18%
- Do not use if the inner side is polluted
- Do not use if the mask is out of shape
- Do not wash and reuse the mask
- In case of pregnant woman, respiratory and cardiovascular disease patient, kids, the old and the infirm, please consult your doctor for advice if breath is not comfortable after using

CERTIFICATIONS

CE CERTIFICATION



EU TYPE EXAMINATION CERTIFICATE

Certificate No: 2163-PPE-1433
Respiratory protective devices, filtering half masks to protect against particles manufactured by
TOPTEC CO., LTD.
140-22, Cheongdam1-ro, Seodong-myeon, Gwangju-si, Gyeonggi-do, Republic of Korea
are tested and evaluated according to

EN 149:2001 + A1:2009 Respiratory Protective Devices - Filtering Half Masks to Protect Against Particles - Requirements, Testing, Marking

Based on the type examination conducted with the evaluation of test reports, technical file according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 5, it is approved that the product meets the requirements of the regulation.

Product Definition

Brand Name: Air Queen Model: Breeze Mask

Filtering half mask

Classification: FFP2 NR

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with:

- Issuing an appropriate EU Declaration of Conformity according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 9
- Ongoing successful performance in fulfilment of the requirements set out in Personal Protective Equipment Regulation (EU) 2016/425 and harmonised standards, ensured by assessments based on Annex 7 (Module C2) or Annex 8 (Module D) of the regulation no later than 1 year from the beginning of serial production

This certificate is initially issued on 11/09/2020 and will be valid for 5 years, if there is no change in the relevant harmonised standard affecting the essential health and safety requirements.



ŞAHİN KACMAZ
UNIVERSAL CERTIFICATION
Director

Yenişehir Bulvarı Katip Sivas Bina No:44/4

Verify Issued Certificate

Certificate Nr

Company Name (Min First 3 Letters)

SEARCH

Please do not request e-mail verification in case you can verify the certificate with scanning QR Codes. Please compare the information results on the query result on the web page and on the certificate on your side when you make query by QRCode or manual.

The certificate holder, certificate number, standard and model name (if exists) must match.

Query Results

Certificate Holder	TOPTEC CO., LTD.
Certificate Nr	2163-PPE-1433
Certificate Type	EN 149+A1:2009 Module B, EU Type Examination Certificate
Model Name	Air Queen / Breeze Mask FFP2 NR
Valid Through	11 / 09 / 2020
Valid Until	10 / 09 / 2025
Issue Date	11 / 09 / 2020
Status	Valid

CERTIFICATIONS

510K FDA CLEARED



2020

CERTIFICATE OF REGISTRATION

This certifies that:
TOPTEC CO., LTD.
 122 Asanvalley-Ro
 Dunpo-Myeon
 Asan-Si Chungcheongnamdo, KR 310409

is registered with the U.S. Food and Drug Administration for FY 2020 pursuant to Title 21, 807 et seq. of the United States Code of Federal Regulations:

Establishment Registration: **3016790437**
 Device Classification Name: **MASK, SURGICAL**
 Product Code: **FXX**
 Regulation Number: **878.4040**

Registrar Corp
 144 Research Drive, Hampton, Virginia, 23666, USA
 Telephone: +1-757-224-0177 • Fax: +1-757-224-0179

Registrar Corp will confirm that such registration remains effective upon request and presentation of this certificate until the end of the year stated above, unless said registration is terminated after issuance of this certificate. Registrar Corp makes no other representations or warranties, nor does this certificate make any representations or warranties to any person or entity other than the named certificate holder, for whose sole benefit it is issued. This certificate does not denote endorsement or approval of the certificate-holder's device or establishment by the U.S. Food and Drug Administration. Registrar Corp assumes no liability to any person or entity in connection with the foregoing.

Pursuant to 21 CFR 807.39, "Registration of a device establishment or assignment of a registration number does not in any way denote approval of the establishment or its products. Any representation that creates an impression of official approval because of registration or possession of a registration number is misleading and constitutes misbranding."

The U.S. Food and Drug Administration does not issue a certificate of registration, nor does the U.S. Food and Drug Administration recognize a certificate of registration. Registrar Corp is not affiliated with the U.S. Food and Drug Administration.

Registrar Corp
 144 Research Drive, Hampton, Virginia, 23666, USA
 Telephone: +1-757-224-0177 • Fax: +1-757-224-0179
 Info@registrarcorp.com • www.registrarcorp.com

David J. Jansz
 Executive Director
 Registrar Corp
 Dated: May 18, 2020

48020-0001 Registrar 5/18

TOPTEC CO., LTD.	KOREA, REPUBLIC OF	3016790437	2020
<ul style="list-style-type: none"> Mask, Surgical - AIR COOL GREEN MASK, AIR QUEEN, AIR QUEEN BREEZE, AIR QUEEN NANO MASK, Nano Mask, PURE MASK, TECHNO WEB, TECHNOWEB 			Manufacturer

Establishment Registration & Device Listing

[FDA Home](#)
[Medical Devices](#)
[Databases](#)

New Search		Back To Search Results
Proprietary Name:	AIR COOL GREEN MASK; AIR QUEEN; AIR QUEEN BREEZE; AIR QUEEN NANO MASK; Nano Mask; PURE MASK; TECHNO WEB; TECHNOWEB	
Classification Name:	MASK, SURGICAL	
Product Code:	FXX	
Device Class:	2	
Regulation Number:	878.4040	
Medical Specialty:	General & Plastic Surgery	
Registered Establishment Name:	TOPTEC CO., LTD.	
Registered Establishment Number:	3016790437	
Premarket Submission Number:	K172500	
Owner/Operator:	TOPTEC Co., LTD.	
Owner/Operator Number:	10072713	
Establishment Operations:	Manufacturer	

HALAL CERTIFICATIONS



korea halal

한국할랄인증서 CERTIFICATE OF KOREA HALAL

인증번호 (Certification No): KHA-20F-00362904

한국할랄인증원(KHA)은 아래 조직에서 생산되는 제품에 대하여 다음 인증규격의 요구사항에 적합함을 인증합니다.
Korea Halal Authority(KHA) certifies that product manufactured by the Below organization complies with the requirements and certification scope of the following certification standards.

회사명(Name of Company): 주식회사 토펙 / TOPTEC CO., Ltd

회사주소(Company Address): 충남 아산시 둔포면 이산밸리로 122
122, A sanvalley-ro, Dunpo-myeon, Asan-si, Chungcheongnam-do, Korea

인증제품명(Name of Products): 나노 마스크(Air Queen)/ Nano Mask(Air Queen)

인증규격(Certification standard): KHAS 29000-할랄 공산용품 일반 규정
KHAS 29000-General Standards For HALAL Industrial Products

이 인증서에 명시되지 않은 다른 제품들은 할랄 제품으로 인정하지 않습니다.
We do not approve of any other products not specified in this certificate as Halal.

인증발행일자(Certification Date of issue): 29.04.2020

인증만료일자(Certification Expiration Date): 28.04.2021

대표이사
Chief Executive Officer
Safiah Weon-Suk Kim

한국할랄인증원
KOREA HALAL AUTHORITY



NELSON IN/EXHALATION TEST REPORT



A Sotera Health company

Sponsor:
 Kyunghun Chung
 Toptec Co., Ltd.
 140-22, Cheongdam-gang 5-ro, Samsung-myeon
 Gumi-si, Gyeongju-gang-do, 30171
 KOREA, REPUBLIC OF

Determination of Inhalation and Exhalation Resistance for Air-Purifying Respirators Final Report

Test Article: Air Queen Breeze Mask
 Study Number: 1205788-001
 Study Received Date: 04 May 2020
 Testing Facility: Nelson Laboratories, LLC
 6200 S. Redwood Rd.
 Salt Lake City, UT 84123 U.S.A.
 Test Procedure(s): Standard Test Protocol (STP) Number: STP0145 Rev 05
 Deviation(s): None

Summary: This procedure was performed to evaluate the differential pressure of compressed air-purifying particulate respirators in accordance with 42 CFR Part 84.100. The air exchange differential or breathability of respirators was measured for inhalation resistance using NIOSH procedure TEB-APR-STP-0007 and exhalation resistance with NIOSH procedure TEB-APR-STP-0003. The differential pressure technique is a simple application of a basic physical principle employing a water manometer differential (upstream and downstream) of the test material, at a constant flow rate.

According to 42 CFR Part 84.84, pretesting must be performed by all applicants, as part of the application process with NIOSH. Results seen below are part of that pretesting and must be submitted to and accepted by NIOSH for respirator approval.

The inhalation resistance criteria as stated in 42 CFR Part 84.100 is an initial inhalation not exceeding 35 mm water column height pressure. The test articles submitted by the sponsor conform to this NIOSH criterion for airflow resistance.

The exhalation resistance criteria as stated in 42 CFR Part 84.100 is an initial exhalation not exceeding 20 mm water column height pressure. The test articles submitted by the sponsor conform to this NIOSH criterion for airflow resistance.

All test method acceptance criteria were met. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 320.



Sean Shepherd electronically approved for
 Study Director

Curtis Gerow

37 May 2020 17:39 (+00:00)
 Study Completion Date-and-Time



A Sotera Health company

Study Number 1205788-001
 Determination of Inhalation and Exhalation Resistance
 for Air-Purifying Respirators Final Report

Results:

Test Article Number	Inhalation Resistance (mm H ₂ O)	Exhalation Resistance (mm H ₂ O)
1	8.8	8.1
2	8.4	8.3
3	8.2	7.9

Test Method Acceptance Criteria: The resistance measurement for the reference plate must be within ± 3 standard deviations of the mean established in the control chart.

Procedure: A complete respirator was mounted to a test fixture comprised of a metal plate with an approximate 3.5 inch diameter hole in the center to allow airflow to reach the mask. The sample holder was assembled by placing a Plexiglas collar around the test fixture and topping with another metal disc with a 3.5 inch opening in the center. The sample holder is held tightly together with clamps and connected to an air source. The manometer is attached to the sample holder by a connection port on the Plexiglas collar.

Before testing, the manometer was zeroed and the back pressure in the sample holder checked and verified to be acceptable. Resistance measurements were taken with a manometer capable of measuring at least 5 inches of water. For inhalation testing, a negative airflow (vacuum) was applied. For exhalation testing, a positive airflow (compressed air) was used. Airflow was passed through the sample holder at approximately 85 \pm 2 liters per minute (L/min).

NELSON NACL AEROSOL TEST REPORT



Sponsor:
Kyunghun Chung
Toptec Co., Ltd.
140-22, Cheongdamgongup 5-ro, Seongsong-myeon
Gumi-si, Gyeongsangbuk-do, 39171
KOREA, REPUBLIC OF

Sodium Chloride (NaCl) Aerosol Test Final Report

Test Article: Air Queen Breeze Mask
Study Number: 1205789-S01
Study Received Date: 04 May 2020
Testing Facility: Nelson Laboratories, LLC
6280 S. Redwood Rd.
Salt Lake City, UT 84123 U.S.A.
Test Procedure(s): Standard Test Protocol (STP) Number: STP0014 Rev 00
Deviation(s): None

Summary: This procedure was performed to evaluate particulate filter penetration as specified in 42 CFR Part 84 and TEB-APR-STP-0050 for requirements on a N95 respirator. Respirators were conditioned then tested for particle penetration against a polydispersed, sodium chloride (NaCl) particulate aerosol. The challenge aerosol was dried, neutralized, and passed through the test article at a concentration not exceeding 200 mg/m³. The initial airflow resistance and particle penetration for each respirator was determined.

According to 42 CFR Part 84.64, pretesting must be performed by all applicants as part of the application process with NIOSH. Results seen below are part of that pretesting and must be submitted to and accepted by NIOSH for respirator approval.

All test method acceptance criteria were met. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.



Study Number 1205789-S01
Sodium Chloride (NaCl) Aerosol Test Final Report

Results: The NIOSH N95 filter efficiency as stated in 42 CFR Part 84.181 is a minimum efficiency for each filter of 95% (±5% penetration). The test articles submitted by the sponsor conform to the NIOSH N95 criteria for filter efficiency.

Test Article Number	Corrected Initial Airflow Resistance (mm H ₂ O)	Maximum Particle Penetration (%)	Filtration Efficiency (%)
1	14.5	1.50	98.41
2	15.8	1.54	98.46
3	15.4	1.88	98.12
4	13.4	1.03	98.07
5	11.5	3.79	96.21
6	12.2	2.98	97.02
7	12.4	3.40	96.00
8	12.6	2.42	97.56
9	12.6	2.06	97.94
10	12.1	3.77	96.23
11	12.0	3.97	96.03
12	13.2	2.04	97.86
13	14.4	2.12	97.88
14	13.6	2.46	97.54
15	2.0	0.306	99.094
16	11.9	2.43	97.57
17	14.4	2.03	97.97
18	14.2	2.47	97.53
19	12.7	2.14	97.86
20	12.0	2.25	97.75

* The final airflow resistance value for each test article was determined by subtracting out the background resistance from the system.

Test Method Acceptance Criteria: The filter tester must pass the "Tester Set Up" procedure. The airflow resistance and particle penetration of the reference material must be within the limits set by the manufacturer.






Robert Dieker electronically approved for
Study Director

Curtis Gerow

02 Jun 2020 15:37 (+00:00)
Study Completion Date and Time

TUV TEST REPORT

DMT - Prüf- und Ergebnisprotokoll			
Prüfung Corona SARS-Cov-2-Virus Pandemie Atemschutzmasken <small>auf Grundlage Prüfgrundritz Rev. 1 vom 26.03.2020 - erstellt von der DLRG Testing and Certification GmbH und gemäß DIN EN ISO 17025 (IAF) der Deutschen gesetzlichen Unfallversicherung</small>			
Angeborenen			
Berichts-Nr.: 8003019226-1-Rev. 01	Prüfdatum: 13.06. - 16.06.2020	Datum Prüfungsanfrage: 08.06.2020	
Prüfung beantragt durch: TUV NORD CERT GmbH Lagerweg 10 45141 Essen	Endkunde: WardWiz Deutschland GmbH Wandlunger Str. 9 41844 Karben	Unterstützt Profer: 19.06.2020 <i>[Signature]</i> (Tilker)	Unterstützt Leiter Prüfamt: 19.06.2020 <i>[Signature]</i> (Schwarz)
Angaben zum Prüfobjekt		 	
Modell/Typbezeichnung (P): MNS Ausführung: Nasenbügel, Ohrschlaufen Hersteller: TOPTEC Co. Ltd. Materialbeschreibung: Air Queen Nanofiber Filter Mask Markt-/Anwendungsbereich: n. A.			
Lieferumfang			
Verpackungsart: Kunststoffbeutel			
Verpackungseinheit: 1 Stück			
Informationsblätter: auf Verpackung gedruckt			
Prüfungen und Test-Ergebnisse			
1. Sichtprüfung		2. Anlegeprüfung	
Kriterium	Bewertung	Kriterium	Bewertung
Verpackungsart:	Kunststoffbeutel	CPA leicht anzulegen	I. O.
Verpackungseinheit:	1 Böck	CPA leicht abzulegen	I. O.
Verpackung staubdicht:	I. O.	Komfort der Kopfbänderung	I. O.
Verpackung staubdicht:	I. O.	Geschäftliche Undichtigkeiten:	nein
Zustand CPA:	I. O.	wahrnehmbare Luftströmung:	nein
Beurteilung:	bestanden	Beurteilung:	bestanden
3. Zustand nach Temperaturveränderung (24 h bei 70°C)		5. Durchlass-Prüfung (Max.-Wert Prüflinge 1-3)	
Zustand CPA:	formstabil, Zustand unbeschadet	Paraffinöl-Konzentration nach Maske - Reingas C ₂ [mg/m ³]:	0,42
Beurteilung:	bestanden	Paraffinöl-Konzentration vor Filter - Rohgas C ₁ [mg/m ³]:	19,85
4. Zustand nach Gebrauchssimulation		Durchlassgrad P [%]:	2,1
Zustand CPA:	formstabil, Zustand unbeschadet	Sollwertabgleich:	Pmax < 6%
Beurteilung:	bestanden	Beurteilung:	bestanden
6. Prüfung CPA mit Ausatemventil			
Prüfung:		Ergebnis:	
Ergebnis Zugversuch am Ventil mit 10 N:		Nur nicht anwendbar	
Funktionsprüfung nach Belastung bei 300 mmHg:		Nur nicht anwendbar	
7. Prüfung des Atemwiderstandes CPA			
Max.-Wert Ausatemwiderstand bei 160 l/min (Sollwert ≤ 3 mbar)		8. Prüfung der Kennzeichnung	
Blickrichtung:	Druckverlust	Kriterium	Prüfung
gerade aus	182 Pa	Herstelleridentifikation:	TOPTEC Co. Ltd.
nach links	-	Materialidentifikation:	MNS
nach rechts	-	Informationsschicht:	auf Verpackung gedruckt
nach oben	-	Aufzeichnung:	auf Verpackung gedruckt
nach unten	-	Hinweise zur Verwendung:	auf Verpackung gedruckt
Max.-Wert Einatemwiderstand bei 95 l/min (Sollwert ≤ 3 mbar)		Beurteilung:	
Blickrichtung:	Druckverlust	bestanden	
gerade aus	165 Pa		
Gesamtbeurteilung der geprüften CPA:		bestanden	
DMT GmbH & Co. KG Am TÜV 1 45127 Essen / Germany TÜV NORD GROUP			

5. Durchlass-Prüfung (Max.-Wert Prüflinge 1-3)	
Paraffinöl-Konzentration nach Maske - Reingas C ₂ [mg/m ³]:	0,42
Paraffinöl-Konzentration vor Filter - Rohgas C ₁ [mg/m ³]:	19,85
Durchlassgrad P [%]:	2,1
Sollwertabgleich:	Pmax < 6%
Beurteilung:	bestanden
Beurteilung:	bestanden

7. Prüfung des Atemwiderstandes CPA		
Max.-Wert Ausatemwiderstand bei 160 l/min (Sollwert ≤ 3 mbar)		
Blickrichtung:	Druckverlust	Beurteilung
gerade aus	182 Pa	bestanden
nach links	-	-
nach rechts	-	-
nach oben	-	-
nach unten	-	-
Max.-Wert Einatemwiderstand bei 95 l/min (Sollwert ≤ 3 mbar)		
Blickrichtung:	Druckverlust	Beurteilung
gerade aus	165 Pa	bestanden
Gesamtbeurteilung der geprüften CPA:		

SGS REPORT

SGS

Test Report No. F650161/LF-CTSAYHA20-11682 Issued Date: 2020. 09. 03 Page 1 of 2

TOPTEC CO., LTD
122 Asanvalley-ro, Dunpo-myeon
Asan-si, Chungcheongnam-do
Korea

The following sample(s) was/were submitted and identified by/on behalf of the client as:-

SGS File No. : AYHA20-11682
Product Name : AIRQUEEN BREEZE MASK
Item No./Part No. : N/A
Received Date : 2020. 08. 31
Test Period : 2020. 08. 31 to 2020. 09. 03
Report Comments : The test part was selected by client's request.
Test Results : For further details, please refer to following page (s)

SGS Korea Co., Ltd.



Tommy Oh / Chemical Lab Mgr

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/sgs-cs/sgs-conditions-of-service> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/sgs-cs/sgs-conditions-of-service>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained herein reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not constitute any legal liability or discharge from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated, the results shown in this test report refer only to the sample(s).

QCP-7081-F01 (01)

SGS Korea Co., Ltd. 305, The Gateway, 76, LG-ro, Daejeon-gu, Arang-m, Gyeonggi-do, Korea 14117
T+82 (0)1 4900 000 F+82 (0)1 4900 004 www.sgs.com/sgs-cs

Member of the SGS Group (Societe Generale de Surveillance)

SGS

Test Report No. F650161/LF-CTSAYHA20-11682 Issued Date: 2020. 09. 03 Page 2 of 2


Sample No. : AYHA20-11682.001
Sample Description : White fabric (outer+inner+filter)
Item No./Part No. : N/A
Materials : N/A

DMAc (N,N-Dimethylacetamide)

Test Item(s)	Unit	Test Method	MDL	Results
DMAc (N,N-Dimethylacetamide)	mg/kg	Ultra sonic extraction with Methylene Chloride / GC-MS	10	N.D.

NOTE: (1) N.D. = Not detected (< MDL)
(2) mg/kg = ppm
(3) MDL = Method Detection Limit
(4) - = No regulation
(5) "" = Qualitative analysis (No Unit)
(6) Negative = Undetectable / Positive = Detectable
(7) The results shown in this test report refer only to the sample(s) tested unless otherwise stated.
This test report is not related to Korea Laboratory Accreditation Scheme.

Picture of Sample as Received :



AYHA20-11682.001

*** End of Report ***

This document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at <http://www.sgs.com/sgs-cs/sgs-conditions-of-service> and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at <http://www.sgs.com/sgs-cs/sgs-conditions-of-service>. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained herein reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not constitute any legal liability or discharge from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated, the results shown in this test report refer only to the sample(s).

QCP-7081-F01 (01)

SGS Korea Co., Ltd. 305, The Gateway, 76, LG-ro, Daejeon-gu, Arang-m, Gyeonggi-do, Korea 14117
T+82 (0)1 4900 000 F+82 (0)1 4900 004 www.sgs.com/sgs-cs

Member of the SGS Group (Societe Generale de Surveillance)

TOPTEC

CERTIFICATIONS

ISO 9001



ISO 14001



OHSAS 18001



ISO 13485



Certification Name	Certification Contents	Certification Date	Remarks
ISO 9001	Quality management system	Jul. 2017	
ISO 14001	Environment management system	Jul. 2017	
OHSAS 18001	Occupational health and safety management system	Jul. 2016	
ISO 13485	Design, Manufacturing and Sales of Masks	Jul. 2020	



Certificate No : RMDM0013

Certificate of Registration
QUALITY MANAGEMENT SYSTEM

This is to certify that
the quality management system of

TOPTEC Co. Ltd.,

at

122, Asanvalley-ro, Dunpo-myeon, Asan-si, Chungcheongnam-do, Korea

Has been found to conform to the Quality Management System Standards:

KS P ISO 13485:2016 / ISO 13485:2016

This Certificate is valid for the following product or service ranges:

Design, Manufacturing and Sales of Masks

Issue Date : Jul. 27. 2020

Certification Date : Jul. 27. 2020

Valid Date : Jul. 26. 2023



Authorized By

Eun-Ju Hwang, President



Quality Management System Certificate

GERMAN CERT

TOPTEC Co., Ltd.

#122 Asan valley-ro, Dunpo-myeon, Asan-si, Chungcheongnam-do, Korea

German Cert Co., Ltd. Hereby certifies that the Quality Management System of the above organization has been evaluated and found to be in line with the requirements of the following standard:

ISO 9001:2015 / KS Q ISO 9001:2015

For the scope of

Design, Development and Manufacture of Automatic Facility and Machine

Certificate Number : **GCQ-202035**

Initial Certification Date : 30 July 2005

Certification Date : 29 July 2020

Expiry Date : 28 July 2023

Issue Date : 24 July 2020

Daek Wook Ki

Scheme Manager



GERMAN CERT Co., Ltd. is accredited by Korea Accreditation Board as a certification body of Quality management system.
(Accreditation number : KAB-QC-08)



GERMAN CERT
QUALITY MANAGEMENT SYSTEM



OH&S Management System Certificate

GERMAN CERT

TOPTEC Co., Ltd.

#122 Asan valley-ro, Dunpo-myeon, Asan-si, Chungcheongnam-do, Korea

German Cert Co., Ltd. Hereby certifies that the Occupational Health & Safety Management System of the above organization has been evaluated and found to be in line with the requirements of the following standard:

ISO 45001:2018 / KS Q ISO 45001:2018

For the scope of

Design, Development and Manufacture of Automatic Facility and Machine

Certificate Number : **GCO-200011**

Initial Certification Date : 29 July 2016

Certification Date : 28 February 2020

Expiry Date : 28 July 2022

Issue Date : 24 July 2020 (Rev. 1-TS)

GERMAN CERT
OCCUPATIONAL HEALTH AND SAFETY MANAGEMENT SYSTEM

Daek Wooki

Scheme Manager



GERMAN CERT Co., Ltd. is accredited by Korea Accreditation Board as a certification body of Occupational Health and Safety management system. (Accreditation number : KAB-OC-32)





GERMAN CERT

OH&S Management System Certificate

TOPTEC Co., Ltd.

#122 Asan valley-ro, Dunpo-myeon, Asan-si, Chungcheongnam-do, Korea

German Cert Co., Ltd. Hereby certifies that the Occupational Health & Safety Management System of the above organization has been evaluated and found to be in line with the requirements of the following standard:

OHSAS 18001:2007

For the scope of

Design, Development, Manufacture and Service of Automatic Facility and Machine

Certificate Number : **KorO-160013**

Initial Certification Date : 29 July 2016

Certification Date : 29 July 2019

Expiry Date : 28 July 2022

Issue Date : 26 July 2019

Daek Wooki

Scheme Manager



GERMAN CERT
OCCUPATIONAL HEALTH AND SAFETY MANAGEMENT SYSTEM



Environmental Management System Certificate

GERMAN CERT

TOPTEC Co., Ltd.

#122 Asan valley-ro, Dunpo-myeon, Asan-si, Chungcheongnam-do, Korea

German Cert Co., Ltd. Hereby certifies that the Environmental Management System of the above organization has been evaluated and found to be in line with the requirements of the following standard:

ISO 14001:2015 / KS I ISO 14001:2015

For the scope of

Design, Development and Manufacture of Automatic Facility and Machine

Certificate Number : **GCE-200769**

Initial Certification Date : 30 July 2005

Certification Date : 29 July 2020

Expiry Date : 28 July 2023

Issue Date : 24 July 2020

Daek Wooki

Scheme Manager



KAB-EC-33

This accreditation mark indicates that the certification body (accreditation number: KAB-EC-33) is accredited by Korea Accreditation Board (KAB) as an environmental management system certification body.



GERMAN CERT
ENVIRONMENTAL MANAGEMENT SYSTEM