



CR3110 Series 100 gloves/box 10 boxes/case

XS: 3115 S: 3116 **M**: 3117 L: 3118

XL: 3119

Contour

Nitrile Powder Free Examination Gloves

- Contours to hands for reduced hand fatigue during extended wear
- Beaded cuff for extra tear resistance
- Fingertip textured in a blue color

Product Information					
Material	Nitrile		Finger Thickness	0.17 mm / 5.2 mil	
Color	Blue		Palm Thickness	0.15 mm / 5.0 mil	
Length (mm)	240		Cuff Thickness	0.10 mm / 4.9 mil	
Texture	Fingertip		Powder Content	Powder-Free	
Sterile	No		Tensile Strength	Unaged: 21 MPa Aged: 16 MPa	
Product Sizing					
Palm Width (mm)	XS (CR3115)	Small (CR3116)	Medium (CR3117)	Large (CR3118)	XL (CR3119)
	70 - 79	80 - 89	90 - 99	100 - 109	110 - 120

Ecological Information

Bio-Compatibility: The chemical formulation of the gloves and surface lubricating materials do not contain any substances normally known to be harmful to the user or to any person with whom the gloves contact.

Other information

To prolong product life, gloves should be stored in a cool, dry location, located away from heat and light.

The data in this Material Safety Data Sheet relates to the specific material designated herein and does not relate to the use in combination with other material or in any process.

No Breakthrough after 240 minutes

No Breakthrough after 240 minutes No Breakthrough after 240 minutes

No Breakthrough after 240 minutes

No Breakthrough after 240 minutes

No Breakthrough after 240 minutes

ASTM D6978 Chemotherapy Drugs Permeation Test

*denotes mandatory chemotherapy drugs

Carmustine* (BCNU) (3.3 mg/ml)

Cisplatin (1.0 mg/ml)

Cyclophosphamide (Cytoxan)* (20.0 mg/ml)

Cytarabine (100 mg/ml)

Dacarbazine (DTIC) (10.0 mg/ml)

Doxorubicin Hydrochloride* (2.0 mg/ml)

Etoposide* (20.0 mg/ml)

Ifosfamide (50.0 mg/ml)

Methotrexate (25.0 mg/ml)

Mitomycin C (0.5 mg/ml)

Mitoxantrone (2.0 mg/ml) Paclitaxel (Taxol)* (6.0 mg/ml)

Thiotepa* (10.0 mg/ml)

Vincristine Sulfate (1.0 mg/ml)

Distributed by:

CRANBERRY USA INC.

Fluorouracil* (50.0 mg/ml)

No Breakthrough after 240 minutes No Breakthrough after 240 minutes

Not Recommended

No Breakthrough after 240 minutes

No Breakthrough after 240 minutes No Breakthrough after 240 minutes

No Breakthrough after 240 minutes

Not Recommended

No Breakthrough after 240 minutes

