

# 3M™ NITRILE GLOVES

AQL 1.5(Compliance Standards : ASTM D6319/EN455 Part 1,2,3&4, FDA 510K)



## 3M Powder-Free NITRILE Medical Gloves MNT760PF



Weight : 4.4g ± 2%



*Medical*



*Veterinary clinic*



*Laboratory*

- Low in Chemical Residuals supported by Electrochemistry Data under IEST-RP C005-2 for cleanroom requirements. Reduced chemical irritation is goal.
- Gloves are easily put on and soft for comfortable wear. high modulus technology makes the gloves odor-free and non-tacky.
- The length of 12-inch gives the user essential protection in the medical use area as a patient examination glove
- The NITRILE glove, is one of the best materials to replace Latex. avoiding the potential for Latex-Sensitive individuals who are known to have allergic reactions and to avoid sensitizing individuals through prolonged use of latex. The Nitrile material is 100% Nitrile. It is a POWDER-FREE Glove.

# 3M™ NITRILE GLOVES

AQL 1.5(Compliance Standards : ASTM D6319/EN455 Part 1,2,3&4, FDA 510K)




## 3M Powder-Free NITRILE Medical Gloves MNT760PF



Weight : 4.4g ± 2%

### 1. Width & Length 12 Inch Glove

## SIZING CHART



Product Sizes	Width (Averages)	Length (Averages)
SMALL	85 +/- 10mm	300mm minimum
MEDIUM	95 +/- 10mm	300mm minimum
LARGE	106 +/- 10 mm	300mm minimum
X-LARGE	116 +/- 10mm	300mm minimum

- Length of all gloves shall be not less than 300 mm for Medium and Large, and not less than 290 mm for Small.

### 2. Nitrile Gloves User Instructions



This is a powder-free glove that does not contain powder.



This product is non-sterile.



There is a fingertip pattern design on the finger part to prevent slipping.



This product is disposable and cannot be reused.



Both hands can be worn regardless of left and right.

### 3. Packaging & Shipping



- Packaging : 100 pieces gloves per box
- Packaging Detail : 100pcs/box, 10box/carton.  
1,035 carton / 20ft container, 2,070 carton / 40ft container

CARTON DIMENSION [ 435 \* 260 \* 250 mm ]